
Sorting single-walled carbon nanotubes by their physical and electronic structure using density differentiation



Mark C. Hersam

Professor of Material Science and Engineering Northwestern University Evanston, IL 60208-3108, USA Tel: 847-491-2696, Fax: 847-491-7820

E-mail: m-hersam@northwestern.edu http://www.hersam-group.northwestern.edu/


3rd NASA-NIST Workshop on Nanotube Measurements

Gaithersburg, Maryland September 26, 2007

Outline

- Motivation and background information
- Density gradient centrifugation of DNA encapsulated SWNTs
- Density gradient centrifugation of surfactant encapsulated SWNTs
- Applications and commercialization






Outline


- Motivation and background information
- Density gradient centrifugation of DNA encapsulated SWNTs
- Density gradient centrifugation of surfactant encapsulated SWNTs
- Applications and commercialization





Carbon Nanotube Polydispersity


- Many chiralities are produced during current methods of SWNT synthesis.
- Chirality (twist and diameter)
 determine electronic and optical
 properties.
- Bulk quantities of carbon nanotubes pure in one chirality would likely enable:
 - > SWNT optical amplifiers
 - > SWNT thin film transistors
 - > SWNT transparent conductors
 - > SWNT multi-analyte biosensors




Conventional Centrifugation of SWNTs

Nano Letters, **3**, 1549 (2003).


- Raw SWNTs are suspended in an amphiphilic surfactant (e.g., sodium dodecyl sulfate, DNA, sodium cholate, etc.) in heavy water by horn ultrasonication.
- Aggregations of nanotubes are removed by centrifugation.





The nanotube optical absorbance and fluorescence spectra become sharper as bundles are removed by centrifugation.

Conventional Centrifugation Removes Bundles but Polydispersity Remains



- ¹ M.J. O'Connell et al., *Science* **297**, 2002.
- ² S. Bachilo et al., *Science* **298**, 2002.
- ³ A. Hagen et al., *Nano Lett.* **3**, 2003.
- ⁴ J-S. Lauret et al., *PRL* **90.** 2003.

Each peak in the PL spectrum represents a different SWNT chirality.



A clear need exists for sorting SWNTs by their physical and electronic structure

Previous Strategies for Sorting SWNTs

- Metallic-semiconducting isolation
 - Dielectrophoresis¹ (electrical)
 - Chemical reactivity^{2,3} (chemical)
 - Anion exchange⁵ (charge)
 - Controlled electrical breakdown⁴ (electrical)
- Diameter sensitivity
 - Anion exchange⁵ (charge)
- ¹ R. Krupke et al. Science **301**, 2003.
- ² M. S. Strano et al. Science **301**, 2003.
- ³ G. G. Samsonidze, Appl. Phys. Lett. 85, 2004.
- ⁴ P.G. Collins, M.S. Arnold, P. Avouris, Science, 2001.
- ⁵ M. Zheng et al. Nature Mater. 2, 2003.

Desired Attributes of a SWNT Purification Process


- Scalable
- Compatible with SWNTs of all lengths and diameters
- Utilizes non-covalent/reversible functionalization
- Iteratively repeatable
- Economical

We found that none of the pre-existing purification strategies sufficiently met all of these criteria.

Outline

- Motivation and background information
- Density gradient centrifugation of DNA encapsulated SWNTs
- Density gradient centrifugation of surfactant encapsulated SWNTs
- Applications and commercialization





SWNT Density as a Function of Diameter


Nano Letters, 5, 713 (2005).

Without external functionalization/hydration:

$$\rho_{NT} := \frac{4 \, \rho_s}{D}$$


Including external functionalization/hydration:

$$\rho_{NT} := \frac{\rho_s \pi D + \rho_{ext} \pi \left(\left(\frac{D}{2} + t \right)^2 - \frac{D^2}{4} \right)}{\pi \left(\frac{D}{2} + t \right)^2}$$



For small diameters, SWNT density was expected to vary by < 1% between adjacent SWNT diameters.

Separation in Density Gradients



Density Gradient Media

Nano Letters, **5**, 713 (2005).

Table 8.1.	Buoyant	density	(g/ml^{-1})	of	macromolecules	in	density	gradient med	lia
------------	---------	---------	---------------	----	----------------	----	---------	--------------	-----

Macromolecule	CsCl	Cs ₂ SO ₄	Metrizamide	Nycodenz [®]	lodixanol
Native DNA	1.71	1.43	1.12	1.13	1.12
Denatured DNA	1.73	1.45	1.14	1.17	1.16
RNA	>1.9	1.64	1.17	1.18	1.18
Polysaccharides	1.62		1.28	1.29	
Proteins	1.3	1.3	1.27	1.27	1.26
Proteoglycans		1.46			

Data from Ford and Rickwood (1983), Rickwood (1992) and Rickwood and Patel (1996).

DNA-Wrapped Carbon Nanotubes

- Single-stranded DNA
 - 1 mg ssDNA: 1 mg NT: 1 mL 0.1 M NaCl -- ultrasonicate
- Sequences 20-60 bases long, rich in G,C, and T work best
 - We use $GT_{(20)}$






Figure and method from M. Zheng et al. Nature Materials 2003.

Imaging DNA-Wrapped SWNTs







• Scanning probe microscopy images are consistent with the proposed helical wrapping of DNA around the SWNT.

Density Gradient Ultracentrifugation (DGU)


- Centrifugation at 174,000 g(64,000 rpm) for 10.5 hours.
- Density gradient formed using iodixanol and heavy water.
- > DNA-wrapped nanotubes are stable in iodixinol.
- Initially, the SWNTs separate by sedimentation velocity.
- * Then, as they approach their isopycnic point, they separate by density.
- * The density gradient becomes steeper with time.






Resulting Separation


- After 10.5 hours, the NTs have separated into three regions (i-iii).
- (i) Colored bands isolated, fully functionalized nanotubes separated by physical and electronic structure (1.11-1.17 g cm⁻³)
- (ii) Grey colorless region bundled or partially functionalized nanotubes / no separation by structure observed (> 1.17 g cm⁻³)
- (iii) Pellet large aggregates and insoluble material



Fractionation





Resulting Fractions



No Length Dependence

Nano Letters, **5**, 713 (2005).






Similar length distributions were observed throughout regions i and ii (the colored bands and the grey region) of the density gradient.

Quantification of Separation by Diameter

Nano Letters, 5, 713 (2005).

λ_{11s} (nm)	Chiralities	Diameters (Å)
929	(9, 1)	7.57
991	(6, 5), (8, 3)	7.57, 7.71
1040	(7, 5)	8.29
1134	(8, 4), (7, 6)	8.40, 8.95
1199	(8, 6)	9.66
1273	(9, 5), (8, 7)	9.76, 10.32


$$\rho_{NT} := \frac{\rho_s \pi D + \rho_{ext} \pi \left(\left(\frac{D}{2} + t \right)^2 - \frac{D^2}{4} \right)}{\pi \left(\frac{D}{2} + t \right)^2}$$



Separation is most efficient for small SWNT diameters

Separation in a Shallow Gradient

Nano Letters, **5**, 713 (2005).



Spectra before (dotted blue line) and after (solid red line) selecting for CoMoCAT-grown SWNTs 7.6 Å in diameter using density gradient centrifugation.

Optical Properties of Chirality-Enriched SWNTs

Collaboration with Tobias Hertel (Vanderbilt Univ.) and Daniel Resasco (Univ. of Oklahoma)

J. Phys. Chem. C 2007, 111, 3831-3835


3831

Pump-Probe Spectroscopy of Exciton Dynamics in (6,5) Carbon Nanotubes

Zipeng Zhu,† Jared Crochet,† Michael S. Arnold,‡ Mark C. Hersam,‡ Hendrik Ulbricht,§ Daniel Resasco, and Tobias Hertel*,†,⊥

Vanderbilt Institute of Nanoscale Science and Engineering, Department of Physics and Astronomy, Vanderbilt University, Nashville, Tennessee, Department of Materials Science and Engineering, Northwestern University, Evanston, Illinois, Institute for Experimental Physics, University of Vienna, Vienna, Austria, and School of Chemical Engineering and Materials Science, University of Oklahoma, Norman, Oklahoma

Received: October 23, 2006; In Final Form: December 19, 2006



Desired Attributes of a SWNT Purification Process


How does density gradient centrifugation of DNA/SWNTs stand up to our criteria?

- Scalable (YES)
- Compatible with SWNTs of all lengths (YES) and diameters (NO)
- Utilizes non-covalent (YES) / reversible (NO) functionalization
- Iteratively repeatable (NO)
- Economical (NO)

Outline


- Motivation and background information
- Density gradient centrifugation of DNA encapsulated SWNTs
- Density gradient centrifugation of surfactant encapsulated SWNTs
- Applications and commercialization






Surfactants

Nature Nanotechnology, **1**, 60 (2006).




(A) Sodium deoxycholate (B) sodium taurodeoxycholate (C) sodium cholate (D) sodium dodecylbenzene sulfonate (SDDBS) (E) sodium dodecyl sulfate (SDS)




Dispersion and density gradient centrifugation of sodium cholate encapsulated SWNTs of various diameters

Sodium Cholate versus DNA

Nature Nanotechnology, 1, 60 (2006).






Sodium cholate


Single-stranded DNA

- \$0.62/g ------ \$2242.80/g (100 g scale) (largest offered, 150 mg scale)
- Reversible encapsulatation ----- Irreversible wrapping?
- Can solubilize nanotubes _ _ _ _ Most efficiently wraps SWNTs of various diameters _ _ _ _ near 0.7-1.0 nm in diameter

Density Gradient Centrifugation with Sodium Cholate



Sodium cholate encapsulated SWNTs after density gradient centrifugation



Initial material and isolated of the purple layer (dominated by the (6, 5) chirality)

Resulting Fractions

Nature Nanotechnology, 1, 60 (2006).




A wide range of optically pure SWNT samples can be produced in one density gradient centrifugation step


(Raw material = HiPCO SWNTs)


Multiple Iterations of Purification


Nature Nanotechnology, 1, 60 (2006).




- Remove a fraction (in which the concentration of the target chirality is maximized)
- Re-run it in a 2nd gradient
- Repeat












> 97% within 0.02 nm of mean diameter




Mean = 0.76 nm







Photoluminescence intensity (A.U.)



Photoluminescence intensity (A.U.)



Diameter Tunability with Co-Surfactants



Purification of Large Diameter SWNTs

Collaboration with Phaedon Avouris (IBM)


- CoMoCAT-grown and HiPCO-grown SWNTs (7-11 Å) are possibly too small for electronics (large Schottky barrier).
- Can we separate tubes in the 11-16 Å range (laserablation)?




Sodium cholate encapsulated SWNTs grown by the laser ablation method after separation in a density gradient




Purification by Diameter and Electronic Type




Electronic Type Sensitivity with Co-Surfactants



Metal versus Semiconductor Separation



Quantification of Electronic Type Separation



Desired Attributes of a SWNT Purification Process


How does density gradient centrifugation of surfactant/SWNTs stand up to our criteria?

- Scalable (YES)
- Compatible with SWNTs of all lengths (YES) and diameters (YES)
- Utilizes non-covalent (YES) / reversible (YES) functionalization
- Iteratively repeatable (YES)
- Economical (YES)


Outline

- Motivation and background information
- Density gradient centrifugation of DNA encapsulated SWNTs
- Density gradient centrifugation of surfactant encapsulated SWNTs
- Applications and commercialization






Possible Applications for Purified SWNTs




- Flexible, transparent logic using films of pure semiconducting SWNTs
- Flexible, transparent conductors using films of pure metallic SWNTs
- Multi-analyte optical / electronic sensors
- Controlled conductivity films for cell growth

Field Effect Transistors from SWNT Films




Metal versus Semiconductor FETs




Semiconducting SWNT films: $\mu > 20$ cm²/V-sec

Transparent Conductors from Purified Metallic SWNT Films



Sheet resistance decreases by ~10x when using metal SWNTs sorted via density gradient centrifugation


Diameter Sorting Enables Metallic SWNT Films with Tunable Optical Absorption




In addition, > 95 % transparency in the IR (2-2.6 μ m) for sheet resistance < 100 ohms/square.

Semi-Transparent Conductive SWNT Stained Glass

- The topmost buoyant bands contain metallic SWNTs sorted by diameter
- Diameter sorting produces dramatic differences in SWNT color



Multi-analyte Optical Biosensors from Purified Semiconducting SWNTs




- Red-shift observed with increasing external dielectric constant.
- Suggests possibility for optically based biosensors using SWNTs.

Heavy water

80

Cardiomyocytes on SWNT Thin Films


Collaboration with Samuel Stupp (NU Medical School)



Cardiomyocytes on carbon nanotube thin film after 7 days

Commercialization: NanoIntegris

http://www.nanointegris.com/




- Start-up company founded in February, 2007.
- NanoIntegris is scaling up the DGU process in an effort to bring high performance SWNTs to the scientific community.
- Product line is called IsoNanotubesTM:

Product	Description
IsoNanotubes-M	Metallic SWNTs
IsoNanotubes-S	Semiconducting SWNTs
IsoNanotubes-D	Uniform Diameter SWNTs

Summary

- Motivation and background information
- Density gradient centrifugation of DNA encapsulated SWNTs
- Density gradient centrifugation of surfactant encapsulated SWNTs
- Applications and commercialization





Acknowledgements

Collaborators

Dr. Phaedon Avouris Prof. Tobias Hertel Prof. Daniel Resasco Prof. Samuel Stupp

Postdocs

Dr. Edward Foley Dr. John Ireland Dr. Alexander Manasson Dr. Gordana Ostojic

Graduate Students

Michael Arnold
Steve Christensen
David Comstock
Norma Cortes
Alex Green
Hunter Karmel
Josh Kellar

Joe Lee
Ben Leever
Matthew Schmitz
Tim Tyler
Michael Walsh
Qing Hua Wang
Nathan Yoder

Undergraduate Students

Jesse Benck
Christopher Chen
Frank Du
James Fakonas
Shengyao Li
Christopher Liman

Funding provided by:

National Science Foundation CAREER, NIRT
NASA Institute for Nanoelectronics and Computing
Army Research Office PECASE, YIP, DURINT, TATRC
Alfred P. Sloan Research Fellowship
ONR Young Investigator Award
Department of Energy IEC
NanoIntegris
Northwestern University NSEC, MRSEC, IBNAM, NCLT