



Overview

 Brief description of light scattering formalism and experiments

 Applications to the solution behavior of responsive block copolymer assemblies

Experimental Setup





Scattering vector:

$$\mathbf{q} = \mathbf{q}_s - \mathbf{q}_i$$

$$q = |\mathbf{q}| = \frac{4\pi n}{\lambda} \sin\left(\frac{\theta}{2}\right)$$

Light is scattered by local fluctuations in the medium


Considerations


- Where do local fluctuations come from?
 - Density
 - Concentration
 - Osmotic pressure
 - Longitudinal (phonon propagation)
 - Orientation
- Timescale dynamic vs. static
- Lengthscale relative to q

Dynamic Light Scattering

- Commonly used for particle sizing in dilute solution
- Key technique for the characterization of smallmolecule and polymeric surfactants, as well as micelles that are responsive to solution conditions
- Measure time dependence of the scattered intensity via an autocorrelation function

Correlation Functions





Berne and Pecora, <u>Dynamic Light Scattering</u>, 1976, Dover.

 For diffusion, the relaxation time is on the order of tenths of milliseconds Intensity vs. time looks like noise.

If the time between samples is on the same time scale of the dynamics, then the intensities are **correlated**.

$$g^{(2)}(\tau) = \frac{\langle I(t+\tau)I(t)\rangle}{\langle I\rangle^2} = 1 + F(A) \left[g^{(1)}(\tau)\right]^2$$

Contains dynamic information!

$$g^{(1)}(\tau) = e^{-\Gamma \tau} \qquad \Gamma = D_M q^2$$


For diffusive motion: Plot Γ vs. q^2 : Slope equal to D_m

Decay function analysis (I)

Method of Cumulants:

Assume narrow, unimodal distribution:

$$g^{(1)}(\tau) = \sum_{i=1}^{n} G(\Gamma_i) \exp(-\Gamma_i \tau_i)$$



$$\ln g^{(1)}(\tau) = \left[\ln g^{(1)}(0)\right] - \langle \Gamma \rangle t + \frac{\mu_2}{2!} t^2 - \frac{\mu_3}{3!} t^3 + \frac{(\mu_4 + 3\mu_2^2)}{4!} t^4$$

$$\langle \Gamma \rangle = \langle D_M \rangle q^2$$

$$\langle D_m \rangle \approx \langle D_o \rangle = \frac{kT}{6\pi\eta \langle R_h \rangle}$$


$$pdi = \frac{\mu_2}{\left\langle \Gamma \right\rangle^2}$$

Normalized Variance

Want pdi < 0.3

Decay Function Analysis (II)

Model the decay as a distribution of relaxation rates



Depolarized DLS

 Measure the intensity of scattered light polarized perpendicularly to the incident light:

Now:

$$g^{(1)}(\tau) = e^{-(D_m q^2 + 6\theta)\tau}$$


and:

$$\theta = \frac{kT}{8\pi\eta R^3}$$
 Is the rotational diffusion coefficient.

Note the stronger R dependence, but spherical molecules have zero depolarized intensity

By incorporating shape dependences into D and θ , you can extract a size and axial ratio of the particle

Relaxation Functions in DLS



Static Light Scattering

- Measure the time averaged scattering intensity
- Osmotic fluctuations come about from solution thermodynamics
- Obtain information about weight averaged molecular weight, size, shape and solution thermodynamics
- Applications: M_w determination of homopolymers, CMC and aggregation number determination of micelles, shape analysis

SLS Formalism

Based on the Zimm equation:


$$\frac{Kc}{R_{\theta}} = \frac{1}{M_{w}P(q)} + 2A_{2}c + \dots$$

$$K = \frac{4\pi^2 n^2 (\partial n / \partial c)_0^2}{\lambda_0^4 N_A}$$

$$P(q) = \left[\frac{3}{(qR)^{3}} \left(\sin qR - qR\cos qR\right)\right]^{2}$$

Form factor for sphere

A₂ is the second osmotic virial coefficient




Small Angle Light Scattering

- SALS is another type of (depolarized) SLS experiment
- Similar formalism to SAXS/SANS
 - Fit intensity data to obtain P(q)
 - Relate to distance distribution function
- Since q range is much larger than SAXS/SANS,
 SALS works best for large (~1-10 μm) particles
- Applications droplet formation, large vesicle characterization (ie: Norman et al.), phase separated polymers

Rheo-SALS and Birefringence


 Incorporating rheology (or flow birefringence) into the SALS experiment:

Bicontinuous Microemulsions



The bright streak is caused by elongated structures aligned in the direction of the flow, above a critical shear stress

Shear-induced vesicle formation from lamellar surfactant



Scattering is a result of concentration fluctuations that moderately couple to the flow

In both experiments, the fluctuations are long ranged (on the order of microns)

^{*} From Krishnan, et al. Phys. Rev. Let. 2001, 098301

^{**} From Schmidt, et al. Rheol. Acta. 1999, 486

Light Scattering in Practice


 The Savin group is primarily interested in the solution characterization of pH and temperatureresponsive block copolymer assemblies

Techniques include:

- DLS to determine R_h
- SLS (Kratky plots) to determine qualitative shape
- Couple with TEM, rheology and other spectroscopic techniques

Applications of Light Scattering to Block Copolymer Assembly


\$1B market for block copolymers




Poly(1,4 butadiene) homopolymer


What happens when the blocks are not thermodynamically compatible?


Self-Assembly






Solution Morphologies of Amphiphilic Block Copolymers






Vesicles (**f** < **0.33**)



Wormlike micelles (0.33 < f < 0.55)




Spherical micelles (f > 0.55)

Increasing fraction of hydrophilic block →

^{*}TEM Images from Zhang and Eisenberg. J. Am. Chem. Soc. 1996, 118, 3168.

^{**}PS PAA assemblies from 2% DMF solution

Phase Behavior of Amphiphilic Block Copolymers




• 'Phase diagram' of 1,4-PB-P(Lys)

Morphology dictated by:

- (1) Core chain stretching
- (2) Corona chain crowding
- (3) Surface tension

Stimuli Responsive Structures

Question: By changing the pH or temperature of a solution, can we disrupt the balance between thermodynamic parameters to induce changes in aggregate structure?



Answer: YES!

From Checot, et al. Angew. Chem. Int. Ed. 2002, 1340.

3 things can happen:


- (1) Aggregate swells stretching of core and/or corona chains
- (2) Destruction surface tension can no longer hold aggregate together
- (3) Morphology shifts crowding leads to dramatic curvature shift... Maybe?

Applications: Delivery, tunable viscosity, 'smart' materials

Current Directions

Target Molecule: PB-P(Lys)


- At high pH, the side chains of poly(L-Lysine) are deprotonated and the polymer assumes an α -helical or a β -sheet conformation.
- At physiological pH, the side chains are protonated, resulting in an expanded coil configuration.
- Expansion of the poly(L-Lysine) α -helices results in swelling of the aggregate.



Advantages Over 'Traditional' Polyelectrolytes (PEs)

- Polypeptide-based block copolymers switch between rod-coil and coil-coil
- With PEs, modest change in solution dimensions from protonation/deprotonation
- With polypeptides, there is chain extension from both breaking the tight α-helical structure as well as directionality of pepide bond – not a true random coil structure

Synthesis of PB-P(Lys)



Ueda, K., Hirao, A. and Nakahama, S. *Macromolecules* **1990**, 23, 939-945. Babin, J. *et al., Faraday Discuss.* **2005**, 128, 179-192.


** Characterize with DLS, TEM


pH Response of Deprotected Samples

- Deprotection schemes:
 - HBr, catalytic hydrogenation
 - HF cleavage
- Efficient deprotection, leaving the unsaturation intact

pH Response for Low P(Lys) Fractions

Compositions:







FTIR of protected (a) and deprotected (b) PB₁₀₇-P(Lys)₂₇

* Expect bilayer morphology


D_m of PB₁₀₇-P(Lys)₂₇



pH Response of PB_m-P(Lys)_n from DLS




Effect of Ionic Strength




As [NaCl] increases, the size of the vesicles decreases as expected


Effect of Temperature




Intermolecular β -sheet Formation?



Parallel β-sheet formation within corona chains would cause a flattening of the interfacial curvature



Viscosity Enhancement at High T



pH Response for Low P(Lys) Fraction


- pH response suggests vesicles are swelling at low pH.
 This response is reversible.
- CD confirms α-helical secondary structure at high pH, random coil at low pH
- Increase in vesicle size ~ 35%
- A decrease in the Debye length results in a decrease in vesicle size at low pH
- PB₆₀-P(Lys)₁₃ has limited solubility at high pH

Temperature Response for Low P(Lys) Fraction


- As the temperature in increased at high pH, the aggregates in solution exhibit a dramatic increase in size and corresponding increase in solution viscosity
- •This is potentially due to interchain β -sheet formation, which flattens the interfacial curvature of the assembly
- This change occurs over a very narrow pH range
- The timescale of the change is on the order of 60 min
- 'Local' irreversibility implies that there is a local minimum in the free energy with the b-sheet

pH Response for Higher Fraction of P(Lys)

Compositions:




pH Response of PB₁₀₇-P(Lys)₂₀₀ from DLS




Forms Spherical Micelles:

Note that hydrodynamic size is comparable to the PB_{107} -P(Lys)₂₇ system although the length of the P(Lys) chain is much longer

pH Response of 'Symmetric' PB_m-P(Lys)_n from DLS




Kratky Plots





Static light scattering suggests a change in the form factor with pH for the PB_{107} -P(Lys)₁₀₀ sample, but not for the PB_{107} -P(Lys)₂₇ sample

Rice to Sphere Transition?




Possible Mechanism:





 PB_{107} -P(Lys)₁₀₀ (scale = 200nm)

pH Response of 'Symmetric' Copolymers



 Position of the peak depends on P(Lys) block length

- pH response suggests rod-like micelles are swelling as pH is reduced
- At low pH, the micelle size peaks, then drops off with decreasing pH
- This drop is proposed to be a result of a morphological shift from rice-shaped micelles to spherical micelles
- This is consistent with that would happen if the effective fraction of the hydrophilic block were to increase

Future Directions

- Dual pH response:
 - Acid and basic chains
 - Triblock copolymers
- New hydrophobic blocks
- Dual temperature-pH responsive materials
- Incorporation of corona crosslinking sites
- Biocompatibility

Summary

- Light scattering is a powerful technique for characterization of polymeric materials in solution
- PB-P(Lys) block copolymers have been synthesized and show a rich diversity in morphology with varying composition
- Secondary structure changes within the P(Lys) corona chains induce response to solution conditions
- For PB-P(Lys), aggregate size and potentially morphology can change with solution pH and temperature

Funding: UVM Startup



CD: Prof. Bruce Armitage (CMU)