



Focused project: High-throughput methods to measure interfacial tension

Combinatorial microfluidic systems will be developed to automatically measure interfacial tension in emulsions. The systems will allow multiple feedstocks to be mixed in a range of compositions.

Today's Outline

- Preview of microfluidic method development in the Processing Characterization group.
Major objective: *Analogue of the 4-roll mill.*
Versatile test platform, including interfacial tension.
- Trial project and results.
- Specialized device for measuring interfacial tension.



Processing Characterization

Broad Goal:

Develop microfluidic tools to measure material (and fluid) properties

Example:


Develop a microfluidic analogue for the 4-roll mill.

- Flow modeling
- Device fabrication (as previously discussed)
- Computer control
- A host of applications




Four-roll mill

4-roll mill



- Adjustable flow type: rotational, mixed and extensional flow. Pure shear is impossible.
- Measure flow-induced properties at long material residence times.
- Applications include flow-induced molecular orientation, particle dispersion, drop deformation and coalescence.

Micro-channel analog




Additional features:

- *Easily scalable and manufacturable.*
- ***3D** control of flow and drop positions is possible.*
- *Convenient for microscopy and other high-resolution probes.*
- *A wider range of flow type, including pure shear, is possible.*




Flow Modeling for device design



- 2D Stokes calculations.
- Rotational, shear, extensional and mixed flows selected by boundary conditions.
- In each case, the flow strength at the stagnation point is
 $\sim 0.2 \times \text{narrow channel wall shear rate}$.




Apparatus








Preliminary results, Cross-channel

Drop deformation





$$D = \frac{a - b}{a + b}$$

$$\tau_d = \eta_{\text{eff}} R / \sigma$$



Control ports


Drop position





Computer control

- Adjust pump speeds to regulate flow strength.
- Image analysis feedback and PID control keeps drop(s) at the stagnation point.
- Record drop deformation and position.





Drop deformation: measure **Relaxation** or **Steady** deformation?


Drop shape relaxation time

$$\tau_d = \eta_{\text{eff}} R / \sigma$$

1 Pa s

100 μm

10 mN/m



0.022 s

***Need fast devices,
by minimizing flow resistance and capacitance,***


***or
Measure quasi-steady properties***





Specialized device for measuring interfacial tension

Measure quasi-steady drop deformation, appropriate for rapid time scales.

- Active and passive mixing of drop and continuous phase.
- Control residence time and deformation rate.



Song, Tice & Ismagilov
Angew. Chem. (2003)





Conclusions

- Microfluidic devices have been specially designed for rapid interfacial tension measurements.
- Peripheral results demonstrate the feasibility of:
 - flow modeling for device design / characterization
 - device fabrication
 - computer control and feedback methods.

Deliverables

- Design strategy, fabrication procedures and development tools (including flow modeling and control routines) will be available to members.