



# Project Overview

---

## Combi Informatics System

Polymer Division: Wenhua Zhang, Mike Fasolka, Amit Sehgal


Alamgir Karim and Eric J. Amis

CTCMS: Andy Roosen, James Warren

# Project Goal – Closing the Loop

## Informatics & feedback for:

- Data management, integration
- Automated data entry, retrieval
- Protocols to standardize data formats
- Link to and refine other Combi processes
- Library of scientific combi-data
- Data visualization



# Progress

---

- Database system selection
- Database server setup
- Security / data protection issues
- Basic framework of the database
- Backbone interface programming
- Instrument programming

# PostgreSQL

---

- **Sophisticated Object-Relational DBMS**
- **Supporting almost all SQL (Structured Query Language), including subselects, transactions, and user-defined types and functions**
- **One of the most advanced open-source database available anywhere**
- **Commercial support available**

**FOR MORE INFO...**


**[WWW.POSTGRESQL.ORG](http://WWW.POSTGRESQL.ORG)**

# Software Technology

---


- **Programming Language**
  - Python (Open source, Object oriented)
- **Secured Shell Connection with the Database server**
- **ZOPE selected as Web server**
  - Base on Python, Open source, Easy to maintain, Database connectivity
- **Integration with current data analysis software (mostly IDL)**
- **Consideration of future interface with other softwares**

# Combi Informatics System



# Basic Structure of Database

- Straight forward structure in consideration of automation & expansion
- Must reflect the frame-work of the Combi lab
- Establish basic protocol of Combi experiment set (All stages of Combi processes)
- Large objects, such as image, video, result chart, original data file storage



**Query**

```
select result.image, result.filename, result.type, rparameter.thickness,
minstrument.minname, minstrument.minmodel, sample.sid from result,
rparameter, minstrument, sample where sample.sid=1 and
result.rid=rparameter.rid and measure.rid=result.rid and
```

Generate Query  
Execute  
Save Query  
Load Query

Result	Measure	Minstrument	Sample	Use
<input checked="" type="checkbox"/> File Name	<input type="checkbox"/> Measuring Date	<input checked="" type="checkbox"/> Instrument Name	<input checked="" type="checkbox"/> Sample Id	<input type="checkbox"/>
<input checked="" type="checkbox"/> Image	<input checked="" type="checkbox"/> Thickness	<input checked="" type="checkbox"/> Model	<input type="checkbox"/> Sample Name	<input type="checkbox"/>
<input checked="" type="checkbox"/> Image Type	<input type="checkbox"/> Temperature	<input type="checkbox"/> Function	<input type="checkbox"/> Component	<input type="checkbox"/>
<input type="checkbox"/> Result	<input type="checkbox"/> X	<input type="checkbox"/> Data Type	<input type="checkbox"/> Type	<input type="checkbox"/>
<input type="checkbox"/> Result Type	<input type="checkbox"/> Y			
<input type="checkbox"/> Original Data	<input type="checkbox"/> Comment			

Set Condition: sample.sid=1

minstrument.minname  
minstrument.minmodel  
sample.sid


Insert  
Delete  
Update  
Refresh  
Next 250  
Save  
Save Info

31	18962	xy	None	0.0	AFM	Dimension	1	
32	18967	Z51.006	None	0.0	AFM	Dimension	1	
33	18972	MP132Z.C	None	0.0	AFM	Dimension	1	
34	18977	x	None	0.0	AFM	Dimension	1	
35	18552	a125epf4	None	60.23185	AFM	Dimension	1	
36	18554	a125epf5	None	60.47415	AFM	Dimension	1	
37	18556	a125epf6	2	56.612875	AFM	Dimension	1	
38	18558	a125epf7	3	53.65325	AFM	Dimension	1	
39	18560	a125epf8	None	53.536825	AFM	Dimension	1	
40	18562	a125epf1	None	52.83155	AFM	Dimension	1	
41	18564	a125epf1	None	55.28465	AFM	Dimension	1	
42	18566	a125epf1	None	52.7116	AFM	Dimension	1	
43	18568	a125epf1	None	None	AFM	Dimension	1	
44	18570	a125epf1	None	None	AFM	Dimension	1	
45	18572	a125epf1	None	74.11975	AFM	Dimension	1	
46	18574	a125epf1	None	72.7472	AFM	Dimension	1	
47	18576	a125epf1	None	73.0979	AFM	Dimension	1	
48	18578	a125epf1	None	72.8991	AFM	Dimension	1	
49	18580	a125ep2f	None	72.28725	AFM	Dimension	1	
50	18582	a125ep2f	None	72.101	AFM	Dimension	1	
51	18584	a125ep2f	None	70.96085	AFM	Dimension	1	
52	18586	a125ep2f	None	70.4291	AFM	Dimension	1	
53	18588	a125ep2f	None	69.8476	AFM	Dimension	1	
54	18982	MP132Z.C	None	0.0	AFM	Dimension	1	

Automatic query generation  
easy to use, no SQL knowledge  
required. Search by range, type,  
user, sample, file, etc.

Result table, image  
and other objects can  
pop-up with a click

**Image**




filename: a125epf14.tif  
itype: None  
thickness: None  
minname: AFM  
minmodel: Dimension  
sid: 1

Image object with its conditions  
shows up on user's desktop  
almost instantaneously



## Database




I

II

III

**Instrument Interface Design Strategy**

# Custom Programming for home made, custom modified instrument



Flow Coater 1 Interface Panel

Change User User Name: Wenhua Sample Name: Sample1

Int Gain: 10 Acceleration: 5 Velocity: 40


Prop Gain: 10 Deceleration: 200 Distance: 35000

Conc(%): 2 Component: 90KPS..2,60KPMMA..{

Start Stop Send Query < 10 >

Comment:

Connecting to Combidatabase ...  
Connected  
New user name: Wenhua



Interferometer Interface Panel

Sample Id: GetInfo CreateId

X Number: 10 Y Number: 10 SpectraName: aaa

X Space: 10 Y Space: 10 ResultFile: bbb

X0: 0 Y0: 0

Start Stop Send Query < 10 >

Comment:

Connecting to Combidatabase ...  
Connected  
New user name: Wenhua


Change User User Name: Wenhua User Dir: c:\tmp

**Custom Program for home made instrument that handle both hardware control and Database interface**

- Uniform interface reduce both programming and training time
- Parameters are recorded and send to the database automatically
- A common query interface that enable inter-instrument communication
- Experiment protocol was enforced and data format standardized
- Experiment procedure are logged for future reference.

# Using Active-X provided by instrument Vender

Similer Interface  
control layout as  
other instrument.



ActiveX control  
provided by vender  
and control camera to  
take photos

Database  
connectivity

Custom features  
are added

- Active X control makes custom instrument control easy to achieve
- add-on features can seamlessly integrated into the program

**AFM-CombiDatabase Interface Panel**

Sample Id:

Method of Collection:

Image Number:  Numbers in X:

X0:  Y0:

X Interval:  Y Interval:

X=-2921, Y=15065 ☒ Send to Database

Thickness:  Temperature:  Time:

Comment:

Experiment start...  
New AFM file generated  
This file is "c:\tmp\Z51.00"

User:


AFM data are parsed and convert to commonly known format

Parameter related to the sample and added on the interface to be send to database

Once experiment start, programming check data folder for new data


For Instrument without an activeX interface, a backdoor method have to be used to get the data and send it to the database.

# Inter-Instrument Communication




**Constantly undated centralized the database  
make data from one instrument immediately  
available for other instrument**

- Stage1: Automated data entry and retrieval




- Stage2: Expert system for each instrument



Online Analysis

- Stage3: Closing the loop of the combinatorial method





# Challenges:

- How the database will response when the amount of data increases,  
how to modify the query process to compensate the effect
- Data standard is needed for better inter-instrument communication and easy real time data analysis
- Real time analysis of multi-dimensional data set
- How reliable the automatic data analysis
- How error in data will be handled
- How the database structure will grow when more function are added in the system
- How it interface with other database system

# Summary

- The system is designed to streamline and complete a combinatorial workflow geared at accelerated materials research.
- Based on open source code and designed for transparency, we intend this system to be a model for parties interested in building LRIS and a test bed for the development of informatics standards.
- To date, our progress has been in building infrastructure, in particular a centralized database system that includes connections to the instrumentation and analysis tools in our combi facilities.
- We anticipate challenges in developing adequate high-throughput data analysis/datamining software for our system,