

Laser Fabrication of Combinatorial Libraries for Polymer and Biotechnology Development

Douglas B. Chrisey

Naval Research Laboratory, Washington, D.C.

Monday, April 26, 2004

NCMC-5: Combinatorial Processing & Characterization



Douglas B. Chrisey

Naval Research Laboratory


Washington, DC 20375

(202) 767-4788; chrisey@nrl.navy.mil

Outline


- Brief Introduction to Combinatorial Processing
- Introduction to MAPLE DW
 - Direct Writing Electronics
 - Direct Writing Biomaterials
 - Rapid Prototyping Biological Systems
- Application of MAPLE DW to the Fabrication of Combinatorial Libraries for Polymer and Biotechnology Development
- Conclusions

Combinatorial Libraries for Polymer and Biotechnology Development



New Materials Needed for Next Generation Applications

Useful Devices/Applications




Improvements in Materials Processing Will Result in Improved Capabilities and New Devices

New Methods of Rapid Material Discovery are Needed

- Conventional Approach to Material Discovery Has Been Edison-ian (Slow and Serendipitous Trial and Error Study)

2-D Parameter Space




Combinatorial Synthesis

- Combinatorial Approaches Have Revolutionized Materials Discovery
 - Automated Parallel Synthesis and High Throughput Screening Systems
 - Very Large Numbers of Compositions Can Be Evaluated Simultaneously
 - Dramatically Reduce the Development Cycle — Speeding Commercialization
- Combinatorial Sample Library Is a 2-Dimensional Scan of Elemental Composition
- Environmentally Sound
- Combinatorial Experiments Conveys the Depth of Material Science In a Single Sample


Combinatorial Processing of Polymer and Biomaterial Libraries

3-D Composition Parameter Space




*Small Amounts of Biomaterials That Are Densely Packed
That Can Be Individually Analyzed!*


Introduction to MAPLE DW



Matrix Assisted Pulsed Laser Evaporation Direct Write (MAPLE DW)




Done Under Ambient Conditions!



‘Ribbons’ in MAPLE DW


The Interaction of the Laser with the Ribbon is the Novelty in MAPLE DW. It is Both a Liability and an Asset.

Liability: Ribbons are Difficult to Fabricate.



Ribbon


Asset: Ribbons Effectively “Quantize” the Material Transferred Making MAPLE DW Coatings Highly Reproducible. Each Laser Pulse Deposits an Identical Mesoscopic “Brick” of Electronic Material.



(Beam Cross-Section X Ribbon Thickness)




Film




Transferred Spot Size vs. Fluence

355 nm, 30 ns FWHM




Solute: BaTiO₃ Nanopowder, 150nm dia.
Matrix: γ-Terpineol + Rheological Modifiers




Transferred Spots Smaller Than “Spot Size” Just-Above-Threshold

Fan Out Network for Contact Redistribution




- MAPLE DW Is a Serial Process, But With High Rep. Rate Lasers (100 kHz), Small Lot Size Manufacturing is Possible
- Can Fabricate Electronics, Power, and Communication and I/O Devices As Well As Viable Cell and Biomaterial Patterns




Direct Writing of Electronics

MAPLE DW of Passive Devices


Spiral Inductor




Interconnects




600 MHz Oscillator




YIG Core Inductor




MAPLE
Direct Write




Fractal Antenna



Transmission Lines




Capacitors



Resistors


Laser Direct-Write and Micromachining of Planar Hydrous RuO_2 Ultra-Capacitor Cells



Volume $\sim 2 \times 10^{-5}$ ml

Mass < 0.1 mg/cap




The Weight of a Human Hair



Actual Machined Ultracapacitor Cells

*Our Technique is Unique Since it Can Direct-Write “Inks”
Containing the Active Material (RuO_x) + The Electrolyte (H_2SO_4)*

Multilayer Capacitors




Prototype Chemiresistive Sensor Subsystem



Chemical


Signal



On Off

Acetone Vapor

Polymer/Carbon



Chemiresistor

*Chemical Signal is Rapidly Adsorbed to Chemoselective Polymer
(Detected as Resistance Change) And More Slowly Desorbed!*

A large school of Yellow Tail Jacks is swimming in clear, blue water above a coral reef. The fish are silvery with a prominent dark lateral line and bright yellow tails. They are moving in various directions, some towards the camera and others away. The coral reef below is dark and textured, with some vertical structures visible. The overall scene is a naturalistic underwater photograph.

Direct Writing of Biomaterials



Yellow Tail Jacks

Diagnostic Sensitivity Requirements

Agent	Infective Dose	Agent	Infective Dose
Anthrax	8,000 to 50,000 spores	Smallpox	10-100 organisms*
Brucellosis	10-100 organisms	VEE	10-100 organisms
Plague	100-500 organisms	Viral Hemorrhagic Fevers	1-10 organisms
Q-fever	1-10 organisms	Botulinum Toxins	~70 ng
Tularemia	10-50 organisms	Staph Enterotoxin B	~30 ng



Biomaterial Sensors Needed for Biological Warfare Agents!

Fabricating a Miniature Dopamine Biosensor Using MAPLE DW




Current Increases With Analyte Concentration

Utilize Engineered Tissue Constructs With Transducers For A Bioelectronic Nose




Engineered Tissue Construct: Olfactory Epithelium



MAPLE DW Approach to Laser Transferring Living Cells and Biomolecules

Transfer Biomaterial


- Proteins
- Antibodies
- Enzymes
- Antibiotics
- DNA
- RNA
- Living Prokaryotic and Eukaryotic Cells



Matrix Material


- Polymers
- Hydrogels
- Sol-Gels
- Proteins
- Aqueous Solution

Side View




Assemble Raw Biomaterial Components ? Living Biological Systems

Picoliter-Scale Dispensing of Biological Solutions




- *Preserves Spot Resolution over Large Transfer Distances (Microwells, Microchannels, etc.)*
- *Demonstrated Volume Transfers From 10 fL to 10 nL/shot*



Laser Spot Size = 125 μm


Protein Microarray Fabrication Using MAPLE DW



- Biotinylated BSA Microarray
- Substrate: Nitrocellulose-Coated Glass Slide
- Spot Size: 50 μm
- Protein Volume Dispensed per Spot: 100 pL
- Fluorescent Tagging Performed with Cy5-Labeled Streptavidin


635 nm Fluorescence Intensity

Single Element Antibody Microarray



- Anti-BSA Microarray
- Substrate: Nitrocellulose-Coated Glass Slide
- Spot Size ~ 150 μm
- Antibody Volume/Spot ~ 1 nL
- *BSA Successfully Identified out of a Mixed Protein Solution (~10 μg/mL)*

30 Different Antibodies Specific to Cell Signaling Pathway Proteins




- ~ 1 nl Antibody Per Spot
- Protein Solution Came From Cells With All Signaling Pathways “Turned-on” So We Should Essentially See All Spots Light up

MAPLE DW of Active Proteins


- Advantages Over Conventional Techniques
 - Works With ANY Starting Material (Viscosity, Powders, Varying Conditions)
 - Reproducibility
 - CAD/CAM
 - Reduced Spot Size $<30 \text{ } \mu\text{m}$
 - Increase Material Utilization/Efficiency by 10^5
 - Works on Planar Substrates and Microwells
 - One System Does Complete Array of Multiple Array Elements

Cell Transfer by MAPLE-DW



Preparation


Culture



Transfer Movie

Transfer

< 1 msec




Analysis

Culture



Single-Cell Resolution and Deposition

Line of Single Human Osteoblasts by Laser Transfer




Do We Genetically Alter the Cells?


- We Use an Arf Excimer Laser at 193 nm = ? (UV) That Is Known to Damage DNA When Directly Exposed to Mammalian Cells.
- Therefore, It Is Possible That Diffuse Laser Light During Our Processing Could Induce Single or Double Strand Break in the Cell's DNA
- We Performed a Neutral Comet Assay on a Group of Transferred P19 Cells and Found No Evidence of Double Strand Breaks in Those Cells
- Alkaline Comet Assay, Which Is More Sensitive to Dna Damage, Also Showed No Damage

Comet Test Results on Laser-Transferred P19 Cells

Day 2



Day 2




- Distributions for 1 hr and 1 Day Post-Transfer Are Nearly Identical To Control Within Statistical Noise

Minimal to No Double Strand Breaks in Laser-Transferred Cells

Rapid Prototyping Living Biological Systems



What if Tissue Could Be Constructed Cell-by-Cell?



- *Is MAPLE DW Suited For Cell-by-Cell Tissue “Construction”?*


- Rapid, Computer-Controlled Placement of Different Materials at 10 to 100 ?m
- Multiple Cell Types
- Molecules Like Growth Factors, Recruitment Factors, Differentiation Chemicals
- Novel Scaffolding Materials (Inorganic/Organic Composites)
- Vascularization (Constructs >1 mm)

Two Cell Type Microarray

Osteoblasts


Cardiocytes

Note the Difference
In Cell Morphology




Mouse Pluripotent Cell Differentiation Post-Laser Transfer

Using DMSO



Muscle Fiber Growth

Using RETINOIC ACID




Observed Dendrite Growth


Transferred Pluripotent Cells Differentiate!

*Next Step to Demonstrate Differentiation
Using Surface Antigen and Cytoskeleton Markers*

Single Shots and Multiple Shots of Rat Cardiac Cells



Single Shot; After 4 Days Culture
Spread Over 700 Microns; 200
Micron Original Spot Size




3 Shots; After 4 Days of
Culture; Cells Did Not Spread;
Bound in Cage/Matrix


*Rat Cardiac Cells Behaved Differently Depending on the
Local Environment!*

Forming Muscle-Like Structures Using MAPLE DW: Lines of Mouse Myoblasts

- 24 hr Incubation Post-Transfer
- Pure Matrigel Substrate



- 48 hr Incubation Post-Transfer
- Pure Matrigel Substrate



Mouse Myoblast Cell Spots Self-Formed Into 3-D Structures!

MAPLE DW of Human Dermal Fibroblasts




Single Cell



2-D Layer

Role of Laser Micro-Patterned Differentially Adherent Scaffolds in Aligning Myoblasts

Randomly Oriented C2C12 Muscle Cells In Culture




Myotube Channels



200 μ m

Agarose
(Non-Adherent)




Matrigel
Matrix
(Adherent)


Solublized Basement
Membrane Preparation:
(Laminin, Collagen IV, Heparan
Sulfate Proteoglycans and Entactin)

Controlled 2-D Adherence Forces Well-Defined 3-D Myoblast Growth?

CAD/CAM Alignment of 2-D Monolayers of Myoblasts to Control Myotube Orientation





Muscle Tissue



*Control of Myotube Size and Formation By Differential Adherence
and Matrigel Channel Width*


Circular Myotube Fabrication Using 2-D Differential Adherence





Controlled Fabrication of Circular Muscles!

Preliminary Tissue Section Transfers By MAPLE DW


- Form Tissue Microarrays by Punching-Out 10 to 1000 μm “Cookies” of Fixed Tissue (Live too?) and Depositing Them Onto Functional Glass Slides
 - Used Gold/Paraffin-Embedded Tissue Layered Structure
 - Used Agarose/Tissue and Agarose/Paraffin-Embedded Tissue Layered Structures



Gold Layer Is Sacrificial



*Agarose Acts to Bind
Tissue “Cookie” Together*




MAPLE DW Versus Other Biomaterial Transfer Techniques


	Laser Capture Microdissection	Microcontact Printing (?CP)	Ink Jet	MAPLE DW
Capable of Patterning Living Cells on Substrates		●		●
Capable of Patterning Active Proteins on Substrates		●	○	●
Capable of Transferring Different Biomaterials Side-by-Side	●		○ <i>*with multiple dispensers</i>	●
Optical Imaging	●			●
Selectively Captures and Immobilizes Cells	●			●
CAD/CAM Compatible			○	●

Engineering Tissue-Based Structures and Devices Cell-By-Cell

- 3-D *Ex Vivo* Tissue Constructs (e.g., Hepatic, Immunological, Neural Living Systems)
- Tissue-Based Sensing of Warfare Agents and Environmental Toxins
- Cell Signaling Platform with Protein Identification Capability
- Cell Separation for Sub-Culturing by Living Microdissector
- Living Microfluidic Devices and Hybrid Biological Motors (Pumps, Valves, Dispenser for Microfluidic Chips)
- Battlefield Repair of Wounded Tissue (*In Situ* Tissue Re-Construction), Computer-Aided Surgeon



Application of MAPLE DW to the Fabrication of Combinatorial Libraries for Polymer and Biotechnology Development



Qualities of MAPLE DW for Combinatorial Processing

- Sub-Picoliter Voxels of Material
- CAD/CAM Process, 100 kHz Deposition Rate, Micron Resolution
- Capability to Work With All Polymers and Biomaterials
- Ability to Deposit Cells One-by-One
- Ability to Deposit 3-D Heterogeneous Structures
 - Engineered Tissue Constructs
 - Tissues Not Found in Nature
 - Living Microdissection
- .

CAD/CAM Tissue Culture Applications

Intracellular Activity:


DNA transcription, protein synthesis, energy metabolism, drug metabolism, cell cycle differentiation, apoptosis

Intracellular Flux:

RNA, hormone receptors, metabolites, calcium, signal transduction, membrane trafficking

Environmental Int.:

Infection, drug action, ligand receptor interactions, cytotoxicity, mutagenesis, carcinogenesis



Cell Products:

Secretion, biotechnology, bioreactor design, product harvesting, downstream processing

Genetics:

Genetic analysis, transfection, infection, transformation, immortalization

Cell-Cell

Interactions:


Morphogenesis, paracrine control, cell proliferation kinetics, metabolic cooperation, cell adhesion and motility, matrix interaction, invasion

Fabricate a Library of Different Polymer Coatings to Test Dissolution Profiles for Controlled Drug Release




**Polymer Library in
Dissolution Environment**


Apply a Specific Biomolecule Drug to a Library of Different Cells



Apply a Library of Different Biomolecule Drugs to Any Number/Type of Cells



Engineered Tissue Construct Surrounded by Antibodies to Cell Signaling Proteins



Combinatorial Libraries by MAPLE DW: Conclusions

- New Era in Laser Processing of Polymers and Biomaterials.
- Parallel Investigation of Various Different Biomaterial/Biomolecule/Cell/Construct/Tissue Combinations.
- MAPLE DW Has Great Potential to Fabricate Unique Polymer and Biomaterial Combinatorial Libraries.