

Multilens Contact Adhesion Tests




Alfred J. Crosby

**Department of Polymer Science & Engineering
University of Massachusetts - Amherst**

NCMC – Adhesion Workshop

October 7, 2002

Outline



- Conventional Adhesion Tests
- Multilens Contact Adhesion Test Theory and Method
- Examples and Implementation
- High-Throughput Analysis
- Practical Points
- Summary


Understanding Polymer Adhesion

Motivation

- “Bolts and screws can be modeled with software... , but glue makers have yet to come up with a predictive model”, **Forbes, 10.29.01**
- Myriad of variables control adhesion
- Existing techniques




Aerospace




Electronic Packaging



Automotive



Biomedical



- Surface Energy
- Molecular Weight
- Time
- Temperature
- Humidity
- Roughness
- Geometry

Peel Tests

Advantages

- Simulates typical use
- Easy sample preparation
- Semi-quantitative results
- Customer-friendly results

Disadvantages

- Results are not absolute
- Sample preparation difficult to standardize
- High statistical populations required
- Stress distribution complicates analysis (backing, adhesive, peel front)



Image from
www.quadgroupinc.com



Axisymmetric Adhesion Tests

Advantages

- Removes backing influence
- Standardizes sample prep
- Bulk and interfacial contributions can be decoupled

Disadvantages

- Absolute analysis is time-consuming
- Consumer knowledge is limited
- Backing properties must be measured separately




Typical Probe-type Adhesion Tests

General Procedure

1. Place adhesive on rigid substrate
2. Position probe above adhesive
3. Move probe into contact with adhesive
4. Hold in contact for arbitrary time
5. Separate probe from adhesive

Typical Measurement Data:



Applied Force
Displacement
(Contact Area)



Spherical Probe



Low Confinement

200 μm





Flat Probe

High Confinement



Spherical Probe

Medium Confinement



Analyzing Probe Tests

Total Dissipation, W_{adh}


- Total area under curve normalized by contact
- Correlates with peel energy

Tack Force, P_{tack}

- Maximum tensile force
- Typically associated with instability initiation

Strain at Failure, δ_{fail}


- Maximum tensile strain



- Quantities are very useful for relative measurements
- Standardization is established, ASTM -
- Contact area images can provide great insight

Theory of Johnson, Kendall, and Roberts (JKR)


General Idea



JKR Equation \longrightarrow
$$a^3 = \frac{9R}{16E} \left[P + 3\pi GR + \sqrt{6GRP + (3\pi GR)^2} \right]$$

Details of Analysis

1. Measure radius of curvature, R
2. Record force, P
3. Measure contact radius, a
4. Correlate a & P data
5. Plot a^3 vs. P and either:
 - Use two parameter fit to determine E and \mathcal{G}
 - or -
 - Measure E independently and determine \mathcal{G}
 - or -
 - Record δ in addition to independently determine E and \mathcal{G}
6. Plot \mathcal{G} vs. da/dt to define material property



Special Cases

Finite-Size Corrections (for $a > h$)



Shull, K.R., et al, *Macromol. Chem. Phys.*, 1998, **199**, 489-511.

Crosby, A.J. et al, *Journal of Applied Physics*, 2001, **88**, 2956-2966.





Viscoelasticity Corrections

Lin, Y.Y., et al, *Journal of Applied Physics*, 1999, **32**, 2250-2260.

Johnson, K.L., ACS publication, 2000.




Outline







- Conventional Adhesion Tests
- Multilens Contact Adhesion Test Theory and Method
- Examples and Implementation
- High-Throughput Analysis
- Practical Points
- Summary

A *Combinatorial* Adhesion Test: MCAT



- Measure a , d
 - Determine G
-
- Possible Variables:
 - Temperature
 - Thickness
 - Strain
 - Surface Energy


Qualitative Analysis




- Contact history images give qualitative information.
- Quickly map the relative adhesion of a surface or library.

How do we calculate G ?

$$G = \frac{3(P' - P)^2}{32\pi E \cdot a^3} \cdot f_P(a, h)$$



$$C = \frac{3}{8Ea} = \frac{d\delta}{dP} = \frac{\delta' - \delta}{P' - P}$$





$$G = \frac{2E(\delta' - \delta)^2}{3\pi a} \cdot f_\delta(a, h)$$

- Minor misalignment tolerable
- Each lens has a unique δ_0
- E must be measured independently or report G/E
- Modify for confinement or viscoelasticity

The Instrument


Time, cost, and precision







Main Components

- Imaging System (i.e. microscope)
- Automated x-y stage
- Displacement control (i.e. actuator or micrometer)
- Displacement sensor
- Computer with DAQ card
- Alignment system (optional)
- Load cell (optional)

The Libraries




Outline



- Conventional Adhesion Tests
- Multilens Contact Adhesion Test Theory and Method
- Examples and Implementation
- High-Throughput Analysis
- Practical Points
- Summary

Contact Experiment




Conditions:

$$dd/dt = 1 \text{ mm/s}$$

$$h_{\text{PS strip}} = 30 \text{ nm}$$

$$\text{Temperature} = 25^{\circ}\text{C}$$

Contact at Elevated Temperatures



Conditions:


$$dd/dt = 1 \text{ mm/s}$$

$$h_{\text{PS strip}} = 30 \text{ nm}$$

Temperature $\sim 80^\circ\text{C}$


Qualitative Analysis

(a)




Room Temperature

(b)




Elevated Temperature

(c)




Simulation

Mapping Multivariable Environments




Quantitative Measurements


PS/PS Contact History



PS/PDMS Contact History







Same Sample, Same Conditions!



Time-Temperature dependence for interfacial strength development indicated by changing weld spot size and thickness dependence.

Outline





- Conventional Adhesion Tests
- Multilens Contact Adhesion Test Theory and Method
- Examples and Implementation
- High-Throughput Analysis
- Practical Points
- Summary

Automated Analysis

Combinatorial Adhesion Analysis

<p>Directory Path for Input Files c:\test\</p> <p>Directory Path for Output Files c:\test_out\</p> <p>Path for Composite Image c:\test_out\</p> <p>Threshold Values</p> <p>Minimum Value: 75</p> <p>Maximum Value: 255</p> <p>0 50 100 150 200 255</p>	<p>Number of Files 820</p> <p>Time of First Image 550.183</p> <p>Pixel Conversion 0.5971</p> <p>E (Pa) 1.00E+6</p> <p>R (m) 3.65E-4</p> <p>Current Image 820</p>
--	--




“Jump” Into Contact

JKR “Zero Force” contact radius

$$a_o = \left(\frac{27\pi R^2 G}{8E} \right)^{1/3}$$


For soft, elastic solids, G is defined by thermodynamic work of adhesion.



$\delta_{\text{critical}} = \delta$ at “jump” into contact



Defined by balance of surface and elastic restoring forces

MultiLens Contact Technique



Advantages of Multilens Contact

- Maximize sensitivity and visualize contact
- Investigate effects of roughness on surface interactions
- Investigate dynamics of surface attraction and separation



Materials

Lens Arrays




- Crosslinked Polydimethylsiloxane (PDMS)
- $E = 1.0 \text{ Mpa}$
- $\gamma \sim 20 \text{ mJ/m}^2$

Substrates


- Bare glass slide (ozone cleaned)
 - $\gamma = 67 \text{ mJ/m}^2$
- n-octyl dimethylchlorosilane coated glass slide
 - $\gamma = 28 \text{ mJ/m}^2$
- Fluorinated glass slide
 - $\gamma = 10 \text{ mJ/m}^2$

Quantifying Jump Distance



Approach
substrate
at 4 nm/s



Dynamics of Contact




Bare Glass



n-octyl
dimethyl




Fluorinated



Time
(s)



Outline




- Conventional Adhesion Tests
- Multilens Contact Adhesion Test Theory and Method
- Examples and Implementation
- High-Throughput Analysis
- Practical Points
- Summary

Practical Points



- Every point does not need to be different.
- Onset of cavitation or fingering is related to strain/stress of “tack force”
- Contact adhesion tests can be used for both weak and strong adhesion
- Control G to measure mechanical properties

Practical Points



- Size of lens and array can be modified
- A row of standard contacts can be incorporated
- Either lenses or substrate can be designed with softer materials



Summary

- Quantitative and qualitative adhesion testing is difficult
- Axisymmetric adhesion tests offer enhanced standardization and information
- Combinatorial approaches can simplify screening processes
- MCAT methodology is general